# **IPS e.max® Abutment Solutions**

## **Cem Kit**

### **Instructions for Use**

Page 2

### **Istruzioni d'uso**

Pagina 41

### **Gebrauchsinformation**

Seite 15

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Pagina 54

### **Mode d'emploi**

Page 28

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Página 67

**CE 0123**

For dental use only!  
Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist  
Made in Liechtenstein  
Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

**ivoclar**  
**vivadent**  
passion vision innovation

## Description

The IPS e.max® Abutment Solutions Cem Kit is a professional kit containing high-quality materials for the cementation of hybrid implant restorations made of IPS e.max lithium disilicate ceramic LS<sub>2</sub> (e.g. IPS e.max Press) or zirconium oxide ceramic.

## Components and their composition

The IPS e.max Abutment Solutions Cem Kit contains the following components:

- **Multilink Implant® – self-curing luting composite with light-curing option**

The monomer matrix is composed of dimethacrylate and HEMA. The inorganic fillers include barium glass, ytterbium trifluoride and spheroid mixed oxide. The particle size is between 0.25 and 3.0 µm. The mean particle size is 0.9 µm. The total volume of inorganic fillers is approx. 40%.

- **Monobond® Plus – primer designed for establishing an adhesive bond between luting composites and all indirect restorative materials**

Alcohol solution of silane methacrylate, phosphoric acid methacrylate and sulphide methacrylate

- **Liquid Strip – oxygen-impermeable glycerine gel to prevent the formation of an oxygen-inhibited layer on composite surfaces**

Glycerine gel

- **Virtual® Extra Light Body Fast Set – polyvinyl siloxane-based impression material (addition-curing) for tacking the ceramic meso-structure on the titanium base as well as for tacking in place the crown during try-in**

Addition-curing silicone containing polyvinyl siloxane, hydrogen siloxane, organic platinum complex, silicate and food colouring

## Indication

Permanent luting of hybrid implant restorations.

Tacking of the ceramic meso-structure on the titanium base as well as tacking in place of the crown during try-in.

## Contraindication

Use is contraindicated

- if the stipulated working technique cannot be applied;
- if a patient is known to be allergic to any of the ingredients of Multilink Implant, Monobond Plus and Virtual Extra Light Body.

## Side effects

- **Multilink Implant, Monobond Plus:**

Systemic side effects are not known to date. In individual cases, allergic reactions to individual components have been reported.

## Interactions

- **Multilink Implant:**

Phenolic substances (e.g. wintergreen oil, eugenol) inhibit polymerization. Consequently, the application of materials containing these substances is to be avoided.

Disinfectants with an oxidative effect (e.g. hydrogen peroxide) may interact with the initiator system, which in turn may impair the curing process.

Therefore, the automix syringe must not be disinfected using oxidative agents. Disinfection may be carried out, for instance, by wiping it with medical alcohol.

- **Monobond Plus:**

If the primer dries up on the surface of the restoration before a strong stream of air has been applied, the surface needs to be rinsed with pure alcohol and blown with air.

**IMPORTANT!** In order to create a strong bond, do not clean the zirconium oxide surfaces with phosphoric acid.

To avoid any contamination, Monobond Plus must not be removed from the bottle directly with a brush.

- **Virtual Extra Light Body Fast Set:**

The setting of polyvinyl siloxanes is inhibited by latex gloves.

The surfaces to which the silicone will be applied (meso-structure, titanium base, etc.) must not be touched with latex gloves.

# e.max Press Abutment Solutions

## Optional: Clinical Try-in

### Temporarily securing the meso-structure on the Ti base

Before the abutment or abutment crown is permanently luted to the Ti base, a clinical try-in may be performed. To facilitate the intraoral handling, the components are temporarily attached to one another with silicone material, e.g. Virtual® Extra Light Body Fast Set.

Please observe the following procedure to temporarily secure the components in place:

- Clean the non-pre-treated Ti base and the pressed object (abutment or abutment crown) with steam and blow dry.
- Place the meso-structure on the Ti base (which is screwed to the model analog) and mark the relative position of the components. This facilitates the achievement of the correct position when the parts are subsequently temporarily assembled.
- Seal the screw channel with a foam pellet.
- Insert the Virtual cartridge into the dispenser and remove the protective cap.
- Screw on the mixing tip and attach the Oral Tip to the mixing tip.
- Apply Virtual Extra Light Body Fast Set both to the Ti base and directly into the meso-structure.
- Insert the Ti base into the meso-structure. Observe the relative position of the objects (rotation lock/markings).
- Hold the parts in the correct relative position for 2:30 minutes until the Virtual Extra Light Body Fast Set has set.
- Carefully remove protruding excess material with a suitable instrument, e.g. a scalpel.

Cleaned, non-pre-treated meso-structures (abutment or abutment crown)

Place the abutment or abutment crown onto the Ti base and mark the relative position.

Seal the screw channel with a foam pellet.

Insert the Virtual cartridge into the dispenser, screw on the mixing tip and attach the Oral Tip.

Apply Virtual Extra Light Body Fast Set both to the Ti base ...

... and directly into the meso-structure (abutment/abutment crown).

Insert the Ti base into the meso-structure. In doing so, observe the relative position of the two components (rotation lock/markings). Hold the components in place for approx. 2:30 minutes until the Virtual Extra Light Body Fast Set material has set.

Carefully remove protruding excess material with a suitable instrument, e.g. a scalpel.

Remove excess Virtual Extra Light Body Fast Set material from the screw channel with an instrument.

## Clinical try-in

### Hybrid abutment

**Important note:** Any intraoral checking of the occlusion/articulation and possible adjustments by grinding may only be performed if the objects have been attached to one another by means of Virtual Extra Light Body Fast Set. During try-in, the Virtual material acts as a buffer, particularly if grinding is necessary, and prevents chipping in the transition area between the hybrid abutment and the crown.

Please observe the following procedure for the clinical try-in:

- Have the clean prepared hybrid abutment (temporarily secured) and the matching clean crown ready at hand.
- Remove the temporary restoration.
- Manually screw in the hybrid abutment with the matching screw.
- Check the geometry of the hybrid abutment (e.g. fit, gingival anaemia) with regard to the gingival margin.
- If required, seal the screw channel on the hybrid abutment with a foam pellet.
- **Tip:** Isolate the inner aspect of the crown with glycerin gel, e.g. Try-in paste, Liquid Strip.
- **Place the crown intraorally onto the hybrid abutment to check and adjust the proximal contacts, if necessary. Note: No occlusal functional checks must be performed at this stage.**
- **For the functional check, the crown has to be secured on the hybrid abutment with Virtual Extra Light Body Fast Set. Do not use Try-in paste for this purpose, as this material is not sufficiently resistant against the compressive forces.**
- Insert the Virtual cartridge into the dispenser and remove the protective cap.
- Screw on the mixing tip and attach the Oral Tip to the mixing tip.
- Apply Virtual Extra Light Body Fast Set to the inner aspect of the crown.
- Use your finger to press the crown to the hybrid abutment until the final position has been achieved. Hold the crown in the final position until the Virtual material has set.
- Remove excess Virtual material.
- Check the occlusion/articulation and make required adjustments with suitable grinding instruments (see separate IPS e.max Recommended grinding instruments for ceramics – use in the dental practice). If adjustments have been made by grinding, conduct another polishing cycle or glaze firing.
- Carefully remove the crown from the hybrid abutment and the hybrid abutment (including Ti base).
- Insert the temporary restoration.

Manually screw in the hybrid abutment with the matching screw. Check the geometry of the hybrid abutment (e.g. fit, gingival anaemia) with regard to the gingival margin.

If required, seal the screw channel of the hybrid abutment with a foam pellet.

**Tip:** Isolate the inner aspect of the crown with glycerine gel.

Place the crown intraorally onto the hybrid abutment to check and adjust the proximal contacts, if necessary. **Note: No occlusal functional checks must be performed at this stage.**

Apply Virtual Extra Light Body Fast Set to the inner aspect of the crown.

Use your finger to press the crown to the hybrid abutment until the final position has been achieved. Hold the crown in the final position until the Virtual material has set.

Remove excess Virtual material.

Check the occlusion/articulation and use suitable grinding instruments to make possibly required adjustments.

Carefully remove the crown from the hybrid abutment and remove the Virtual Extra Light Body Fast Set material.

Unscrew the hybrid abutment.

## Hybrid abutment crown

Please observe the following procedure for the clinical try-in:

- Have the cleaned hybrid abutment crown (temporarily secured with Virtual Extra Light Body Fast Set) at hand.
- Remove the temporary restoration.
- **Place the hybrid abutment crown intraorally onto the implant to check and possibly adjust the proximal contacts. Note: No occlusal functional checks must be performed at this stage.**
- Manually screw in the hybrid abutment crown with the matching screw.
- Check the geometry of the hybrid abutment crown (e.g. fit, gingival anaemia) with regard to the gingiva.
- Check the occlusion/articulation and make possibly required adjustments with suitable grinding instruments (see separate IPS e.max Recommended grinding instruments for ceramics – use in the dental practice). If adjustments have been made by grinding, conduct another polishing cycle or glaze firing.
- Carefully remove the hybrid abutment crown (including Ti base).
- Rinse the implant site e.g. with Cervitec® Liquid (antibacterial mouth wash with chlorhexidine) to clean and disinfect it.
- Insert the temporary restoration.

Place the hybrid abutment crown intraorally onto the implant to check and possibly adjust the proximal contacts. **Note: No occlusal functional checks must be performed at this stage.**

Manually screw in the hybrid abutment crown with the matching screw.

Check the geometry of the hybrid abutment crown (e.g. fit, gingival anaemia) with regard to the gingiva.

Check the occlusion/articulation and use suitable grinding instruments to make possibly required adjustments.

Carefully remove the hybrid abutment crown (including Ti base).

# IPS e.max<sup>®</sup> Press Abutment Solutions

## Permanent Cementation Ti base / Meso-structures

Careful preparation of the bonding surface is a prerequisite for an optimum adhesive cementation of the Ti base and the pressed object. The following paragraphs outline the required procedures. The procedure is the same for hybrid abutments and hybrid abutment crowns.

### IPS e.max Press Abutment Solutions (LS<sub>2</sub>)

	Abutment, abutment crown (LS <sub>2</sub> )	Ti base
Blasting	–	The bonding area with Al <sub>2</sub> O <sub>3</sub> at low pressure
Etching	The bonding area with IPS <sup>®</sup> Ceramic Etching Gel for 20 s	–
Conditioning	The bonding area with Monobond Plus for 60 s	
Adhesive cementation	Multilink Implant MO 0	
Covering the cementation joint	Glycerine gel, e.g. Liquid Strip	
Curing	7-minute curing (optional in a light-curing device)	
Polishing the cementation joint	Customary polishers for ceramic/resin materials	

### Abutments made of zirconium oxide (ZrO<sub>2</sub>)

	Meso-structure (ZrO <sub>2</sub> )	Ti base
Blasting	The bonding area with Al <sub>2</sub> O <sub>3</sub> at low pressure	The bonding area with Al <sub>2</sub> O <sub>3</sub> at low pressure
Conditioning	The bonding area with Monobond Plus for 60 s	
Adhesive cementation	Multilink Implant MO 0	
Covering the cementation joint	Glycerine gel, e.g. Liquid Strip	
Curing	7-minute curing (optional in a light-curing device)	
Polishing the cementation joint	Customary polishers for ceramic/resin materials	

### Pre-treatment of the Ti base

To prepare the Ti base for cementation with the meso-structure, please observe the following procedure:

- Observe processing instructions of the manufacturer of the Ti base.
- Clean the Ti base in an ultrasonic cleaner and blow dry or use a steam cleaner.
- Screw Ti base onto a model analog.

- Place the meso-structure on the Ti base and mark the relative position of the components. This facilitates locating the correct position when the parts are assembled at a later stage.
- The emergence profile of the Ti base must not be blasted or modified in any way. To protect the emergence profile, hard modelling wax is applied, as this material can be easily removed later on.
- Seal the screw channel with wax.
- Carefully blast the bonding area with  $\text{Al}_2\text{O}_3$  (50–100  $\mu\text{m}$ ) at low pressure until an even mat surface has been achieved.
- Clean with an instrument and steam cleaner. Make sure that any wax residue is carefully removed.
- After cleaning, any contamination of the bonding surface must be prevented, since contaminations negatively influence the bond.
- Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.
- Seal the screw channel with a foam pellet or wax. Make sure that the bonding surface is not contaminated.

Screw Ti base onto a model analog. Mark the relative position to the meso-structure.

Apply wax to protect the emergence profile. In addition, seal also the screw channel with wax.

Carefully blast the bonding area with  $\text{Al}_2\text{O}_3$  (50–100  $\mu\text{m}$ ) at low pressure until an even mat surface has been achieved.

Clean with an instrument and steam cleaner. Make sure that any wax residue is carefully removed.

Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.


Seal the screw channel with a foam pellet or wax. Make sure that the bonding surface is not contaminated.

## Preparing the pressed object made of IPS e.max lithium disilicate glass-ceramic (LS<sub>2</sub>)

To prepare the pressed object for cementation onto the Ti base, please observe the following procedure:

- Do **not** blast the IPS e.max Press object to prepare it for cementation.
- Clean the pressed object in an ultrasonic cleaner and blow dry or use a steam cleaner.
- After cleaning, any contamination of the bonding surface must be prevented, since contaminations negatively influence the bond.
- To protect the outer surfaces or the glazed areas, wax may be applied.
- Etch the bonding surface with 5% hydrofluoric acid gel (IPS Ceramic Etching Gel) for 20 seconds.
- Subsequently, thoroughly rinse the bonding surface under running water and dry with oil-free air.
- Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.

Do **not** blast IPS e.max Press object

Etch with IPS Ceramic Etching Gel for 20 seconds.

Allow Monobond Plus to react for 60 seconds and dry with air.

## Preparing the meso-structure made of zirconium oxide (ZrO<sub>2</sub>)

To prepare the meso-structure for cementation onto the Ti base, please observe the following procedure:

- Secure the meso-structure on a holder before blasting to facilitate handling (e.g. Microbrush).
- To protect the outer surfaces of the meso-structure, apply wax.
- To better control the blasting procedure, colour the bonding surface of the meso-structure with a permanent marker.
- Blast the bonding surface of the meso-structure with Al<sub>2</sub>O<sub>3</sub> 100 µm at max. 1 bar / 29 psi pressure and at a distance of max. 1-2 cm to the blasting nozzle.
- Clean the meso-structure with the steam cleaner or in an ultrasonic bath (ethanol bath).
- After cleaning, any contamination of the bonding surface must be prevented, as contaminations negatively influence the bond.
- Apply Monobond Plus on the cleaned bonding surface and allow to react for 60 seconds. After the reaction time, dry the remaining residue with water- and oil-free air.

Apply wax to protect the outer surfaces.

Blast the bonding surface.

Allow Monobond Plus to react for 60 seconds and blow dry.

## Cementation with Multilink® Implant

For an optimum bond between the meso-structure and the Ti base, use the self-curing luting composite Multilink® Implant with light-curing option. Read the respective Instructions for Use for more detailed information.

Please observe the following procedure for cementation:

- Keep the cleaned and conditioned components that are to be luted (meso-structure, Ti base) at hand.
- **The subsequent cementation procedure must be carried out quickly and without interruption. The working time of Multilink Implant is 90 (± 15) seconds at 23 °C (± 1 °C) or 73 °F (± 1.8 °F).**
- As a general rule, attach a new mixing tip to the Multilink Implant syringe prior to each use.
- Apply Multilink Implant directly from the mixing tip in a thin layer on the bonding surface of the **Ti base and the bonding surface of the meso-structure.**
- Leave the mixing tip attached to the Multilink Implant syringe until the next application. As the material will cure in the mixing tip, it will serve as a seal.
- Position the meso-structure above the Ti base in such a way that the position markings are aligned.
- Use even and low pressure to join the parts and check the correct relative position of the components when they are in their final position (transition Ti base/meso-structure).
- Subsequently, tightly press the components together for 5 seconds.
- Carefully remove excess in the screw cavity, e.g. with Microbrush or brush, using a rotary movement.
- Remove excess at the transition to the Ti base carefully in its ductile state, e.g. with a foam pellet, while applying slight pressure to hold the components in place.
- Apply a glycerine gel (e.g. Liquid Strip) to the cementation joint to prevent the formation of an inhibition layer.
- After that, completely polymerize the luting composite for 7 minutes in a light-polymerization device.
- **Important: Do not move the objects until the Multilink Implant material has completely cured and hold them in place without allowing for any motion, e.g. with diamond-coated tweezers.**
- After completion of the polymerization, rinse off the glycerine gel with water.
- **Smooth out and polish the cementation joint with rubber polishers.**
- If there are any cement residues in the screw channel, remove them using suitable rotary instruments.
- Clean with a steam cleaner.

Keep the cleaned and conditioned components that are to be luted at hand.

Attach a new mixing tip to the Multilink Implant syringe prior to each use.

Apply Multilink Implant directly from the mixing tip in a thin layer on the bonding surface of the Ti base.

Apply Multilink Implant directly from the mixing tip in a thin layer on the bonding surface of the meso-structure.

Position the meso-structure above the Ti base in such a way that the position markings are aligned.

Join the components using even and light pressure. Subsequently, tightly press the components together for 5 seconds.

Carefully remove excess in the screw cavity, e.g. with a Microbrush or a brush, using a rotary movement.

Remove excess carefully in its ductile state, e.g. with a foam pellet, while applying slight pressure to hold the components in place.

Apply a glycerine gel (e.g. Liquid Strip) to the cementation joint to prevent the formation of an inhibition layer.

Polymerize the luting composite for 7 minutes (optionally in a light-curing unit). **Important: Do not move the objects until the material has completely cured and hold them in place without allowing for any motion.**

After completion of the polymerization, rinse off the glycerine gel with water.

Smooth out and polish the cementation joint with rubber polishers.

If there are any cement residues in the screw channel, remove them using suitable rotary instruments. Do not damage the Ti base.

Completed hybrid abutment and hybrid abutment crown after cementation.

### Additional notes

- Multilink Implant should be processed at room temperature. Cool temperatures may render the material difficult to dispense and mix and prolong working and curing times.
- Apply Liquid Strip cautiously to avoid blending it with or displacing the composite. Leave Liquid Strip on the cement joint until polymerization is complete.
- Polyvinyl siloxanes such as Virtual Extra Light Body Fast Set are chemically resistant. Unpolymerized material may stain clothing.

### Warning

- Prevent any contact of uncured Multilink Implant and Monobond Plus with skin / mucous membrane and eyes. When unpolymerized, these products may cause slight irritation and lead to a sensitization against methacrylates. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Monobond Plus:
  - Highly flammable
  - Irritating
  - Product contains ethanol!
  - Keep away from sources of ignition.
- Virtual Extra Light Body Fast Set:  
If uncured material comes into contact with the eyes, rinse immediately with plenty of water. If irritation persists, consult a physician. If the material comes into contact with the skin, rinse with water and soap.

### Shelf life and storage

- Do not use the individual components of the IPS e.max Abutment Solutions Cem Kit after the indicated date of expiration.
- Storage temperature of Monobond Plus, Liquid Strip and Virtual Extra Light Body Fast Set:  
2–28 °C/36–82 °F
- Storage temperature of Multilink Implant:  
2–8 °C/36–46 °F
- Shelf life: see imprint on syringe/bottle/cartridge
- Monobond Plus and Liquid Strip:  
Close bottle/syringe immediately after use.
- Multilink Implant and Virtual Extra Light Body Fast Set:  
Leave mixing tip on the syringe/cartridge after use. The material will polymerize within the mixing tip, functioning as a natural cap.

**Keep out of the reach of children!**

**For use in dentistry only.**

**Date information prepared:** 01/2012, Rev. 1

### Manufacturer:

Ivoclar Vivadent AG, FL-9494 Schaan / Liechtenstein

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the products for their suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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## **Ivoclar Vivadent Limited**
## **Ivoclar Vivadent, Inc.**